


General Cytology Preservative



- The entire sample is transferred to the preservative vial
- Less flammable, non-toxic and compatible with shipping standards
- Reduced plastic disposable waste
- Morphology consistent with conventional cytology
- Compatible with most HPV/STD tests



CS-F50



- Semi-automated process (50 slides per hour)
- Dual filter technology to minimize the effects of blood, mucus and inflammatory cells

CELLSOLUTIONS GmbH
Halbinselstr. 37
88142 Wasserburg, Germany
Office: +49 8382 942-9901
Service: +49 157 80804380
Email: cseurope@cellsols.com
Website: cellsols.com



BestCyte



- Fully automated digital cytology slide scanner with throughput of 1000 thin-layer slide preparations per day
- Reduces screening time by condensing an entire slide presentation into a few galleries of digitized cells for diagnostic review
- Entire cell deposit is digitized and available for review if needed
- Networked storage for virtual slides
- Web-based software allows access to case review from multiple locations



CS-120

- Positive sample identification
- Automatically measures cellularity
- Output to compatible staining cassettes
- Fully automatic operation after a cell pellet has been centrifuged
- CS30 processes 30 slides/hour
- CS120 processes 120 slides/hour
- Reduced consumables and waste
- Consistent and optimal cellularity



CS-30